

SARS-CoV-2 Antigen Rapid Test (Swab) Package Insert

REF ISCO-ACO502 English

In vitro diagnostic rapid test for qualitative detection of Nucleocapsid Protein antigens of SARS-CoV-9

INTRODUCTION

The novel coronaviruses belong to the β genus. COVID-19 is an acute respiratory infectious disease. People are generally susceptible. Currently, the patients infected by the novel coronavirus are the main source of infection; asymptomatic infected people can also be an infectious source. Based on the current epidemiological investigation, the incubation period is 1 to 14 days, mostly 3 to 7 days. The main manifestations include fever, fatigue and dry cough. Nasal congestion, runny nose, sore throat, myalgia and diarrhea are found in a few cases.

INTENDED USE

The ACROTM SARS-CoV-2 Antigen Rapid Test (Swab) is a *lateral flow* immunoassay for the qualitative detection of Nucleocapsid protein antigens of SARS-CoV-2 in swab specimens from individuals suspected of COVID-19 by their healthcare provider.

The SARS-CoV-2 antigen is able to be detected in the upper respiratory tract during the acute phase of infection. The ACROTM SARS-CoV-2 Antigen Rapid Test (Swab) tests for the SARS-CoV-2 nucleocapsid protein antigen in upper respiratory specimens. While a positive result indicates the existence of the viral antigen, further evaluation of patient history and diagnostics are necessary to confirm these results. A patient may experience concurrent infection with other viruses or a bacterial infection which are not ruled out by a positive SARS-CoV-2 antigen result. Detection of the antigen does not surely certify it as the cause of the disease.

Negative results should be treated as presumptive and confirmation with a molecular assay may be necessary for patient management. Negative results do not rule out SARS-CoV-2 infection and should not be used as the sole basis for treatment or patient management decisions, including infection control decisions. Negative results should be considered in the context of a patient's recent exposures, history and the presence of clinical signs and symptoms consistent with COVID-19.

PRINCIPLE

The ACROTM SARS-CoV-2 Antigen Rapid Test is a lateral flow immunochromatographic assay for the detection of extracted nucleocapsid protein antigens specific to SARS-CoV-2 in swab specimens.

The swab specimens require a sample preparation step. The extracted specimen is added to the sample well of the test cassette to initiate the test. When the specimen migrates in the test strip, SARS-CoV-2 antigens (if present) bind to anti-SARS-CoV-2 nucleocapsid protein conjugated to the indicator and capture particles in the conjugate pad, forming an immune complex. The complex is then captured by the test line on the nitrocellulose membrane as it migrates through the strip.

The presence of two colored lines in the control line region "C" and test line region "T" indicates SARS-CoV-2 Antigen positive. The presence of one colored lines in the control line region "C" indicates SARS-CoV-2 Antigen negative. No appearance of a colored line in the control region "C" indicates an invalid test.

KIT COMPONENTS

•Test cassettes •Sterile swabs •Extraction buffer •Extraction tubes and tips (Optional)
•Workstation •Procedure card •Package insert

STORAGE AND STABILITY

Store the kit at 2-30°C. **DO NOT FREEZE**. The kit is stable through the expiration date printed on the packaging. The test must remain in the sealed pouch until use.

Kit components must be at room temperature (15–30 °C) when used for testing.

WARNINGS AND PRECAUTIONS

- For in vitro diagnostic use only.
- 2. For professional use only.
- This package insert must be read completely before performing the test. Failure to follow directions in package insert may yield inaccurate test results.
- 4. Do not use kit beyond its expiration date.
- Leave test card sealed in its foil pouch until just before use. Do not use if pouch is damaged or open.
- Do not eat, drink or smoke in the area where the specimens or kits are handled.
- Handle all specimens as if they contain infectious agents. Observe established precautions against
 microbiological hazards throughout in the collection, handling, storage and disposal of patient
 samples and used kit contents.
- Wear appropriate personal protection equipment and gloves when running each test and handling patient specimens. Change gloves between handling of specimens suspected of COVID-19.
- 9. Wash hands thoroughly after handling.
- 10. Proper sample collection, storage and transport are essential for correct results.
- 11. Do not reuse the used test cassette.
- 12. Inadequate or inappropriate sample collection, storage, and transport may yield false test results.
- 13. Do not store or test specimens in viral transport media, as it may result in false positive or false negative results.
- $14. \ All \ components \ of this \ kit \ should \ be \ discarded \ as \ biohazard \ waste \ according \ to \ local \ regulations.$

QUALITY CONTROL

Procedural Controls

This test includes a positive procedural control within the test strip and is found at the "C" Control area. A visible line at the marked "C" area indicates the specimen volume has correctly traveled up the test strip and the technician adequately handled the test. A clear or light pink background elsewhere on the test strip should not contain any dark color, which shows that this negative procedural control confirms an adequately-working test.

External Quality Control

Good Laboratory Practice (GLP) compliance dictates the use of external positive and negative controls. ¹ It is recommended to use these controls; however, they are not currently included with this test kit.

SPECIMEN COLLECTION AND HANDLING



Nasal Swab Specimen Collection:

Insert a sterilized swab into a nostril. Do not insert it more than half an inch. Slowly twist the swab, rubbing it along the insides of the nostril for 15 seconds. Using the same swab repeat the collection procedure in another nostril. Withdraw the swab; avoid excess volume and high-viscous nasal discharge.

Caution: If the swab stick breaks during specimen collection, repeat specimen collection with a new

SPECIMEN TRANSPORT AND STORAGE

For best performance, Swab specimens should be tested as soon as possible after collection.

If immediate testing is not possible, and to maintain best performance and avoid possible contamination, it is highly recommended the nasal swab is placed in a clean, unused plastic tube labeled with patient information, preserving sample integrity, and capped tightly at 2-8°C for up to 24 hours prior to testing.

TEST PROCEDURE

All test components should be balanced at room temperature before use.

Specimen Preparation


Only extraction buffer and tubes provided in the kit are to be used for swab specimen preparation.

- Place the collected swab specimen into the extraction tube with extraction buffer. Press the swab against the tube and rotate the swab for 10-15 seconds. (See procedure card for detailed information on specimen Preparation.)
- Remove the swab while squeezing the swab head against the inside of the extraction tube to express as much liquid as possible from the swab. Discard the swab.

*NOTE: Extracted specimen solution is stable for 2 hours at room temperature or 24 hours at 2-8°C.

Test Reaction

- Remove the test cassette from the sealed foil pouch and use it within one hour. Best results will be
 obtained if the test is performed immediately after opening the foil pouch.
- Invert the specimen extraction tube, add 3 drops of extracted specimen (approx.75-100µl) to the specimen well(S) and then start the timer.
- Read the result in the test window 15 minutes after sample application. Results should not be read after 20 minutes.



TEST INTERPRETATION

NEGATIVE: A

A negative specimen will give a single colored line in the control region (C) in the test window, indicating a negative result. This control line means that the detection part of the test was done correctly, but no SARS-CoV-2 Antigen was detected.

POSITIVE:

A positive specimen will give two colored lines. This means that SARS-CoV-2 Antigen was detected. Specimens with low levels of antigen may give a faint Line. Any visible colored line is considered a positive.

INVALID:

If no Control lines are seen, the test is invalid. Insufficient specimen volume or incorrect procedural techniques are the most likely reasons for control line failure. Review the procedure and repeat the test with a new test. If the problem persists, discontinue using the test kit immediately and contact your local distributor.

TEST LIMITATIONS

- The test is intended for direct swab specimens only. Viral Transport Media (VTM) should not lused with this test as it may cause false results.
- 2. When testing for detection of SARS-CoV-2 Nucleocapsid Protein Antigen, it is imperative to carefully follow the provided procedure, results interpretation, and specimen collection for a valid test result. Deviating from these instructions may lead to erroneous results.
- The performance of the ACROTM SARS-CoV-2 Antigen Rapid Test (Swab) was evaluated using the procedures provided in this product insert only. Modifications to these procedures may alter the performance of the test.
- Neither the quantitative value nor the rate of increase in the concentration of SARS-CoV-2 antigens can be determined by this qualitative test.
- 5. The ACROTM SARS-CoV-2 Antigen Rapid Test (Swab) will only indicate the presence of SARS-CoV-2 Antigens in the specimen and should not be used as the sole criteria for the diagnosis of SARS-CoV-2 infections.
- 6. If the test result is negative or non-reactive and clinical symptoms persist, it is recommended to re-sample the patient and test again or test with a molecular diagnostic device to rule out infection in these individuals.
- A negative test result may occur if the level of antigen in a sample is below the detection limit of the test or if the sample was collected improperly.
- The results obtained with the test should be considered with other clinical findings from other laboratory tests and evaluations.
- Negative results do not rule out SARS-CoV-2 infection, particularly in those who have been in contact with the virus. Follow-up testing with a molecular diagnostic should be considered to

rule out infection in these individuals.

10. Positive results of SARS-CoV-2 may be due to infection with non-SARS-CoV-2 coronavirus strains or other interference factors.

PERFORMANCE CHARACTERISTICS

1. Clinical performance

Clinical performance characteristic of the ACROTM SARS-CoV-2 Antigen Rapid Test (Swab) was evaluated with 421 swab specimens. RT-PCR is used as the reference method.

The data was summarized in the table below:

Nasal swab specimen

ACROTM SARS-CoV-2		R'	Total		
Antigen Rapid Test		Positive	Negative	Total	
SARS-CoV-2	Positive	105	2	107	
Antigen	Negative	6	308	314	
Total		111	310	421	
Relative Sensitivity		105/111=94.6% (95%CI*: 88.6%~98.0%)			
Relative Specificity		308/310>99.4% (95%CI*: 97.7%~99.9%)			
Accuracy		413/421=98.1% (95%CI*: 96.3%~99.2%)			

^{*95%} Confidence Intervals

2. Detection limit

 $ACRO^{TM}\,SARS\text{-}CoV\text{-}2\,Antigen\,Rapid\,Test\,was\,confirmed\,to\,detect\,100\text{TCID}_{50}/\text{ml}\,of\,SARS\text{-}CoV\text{-}2.$

3. Cross Reactivity (Analytical Specificity) and Microbial Interference

The following viral strains were tested at concentrations in following table and all found to be negative when tested with the ACROTM SARS-CoV-2 Antigen Rapid Test (Swab):

	Test items	Test Concentration
	Respiratory syncytial virus	8.89 x 10 ⁴ TCID ₅₀ /ml
	Parainfluenza virus 3	1.58 x 10 ⁸ TCID ₅₀ /ml
	Parainfluenza virus 2	1.58 x 10 ⁷ TCID ₅₀ /ml
	Mumps	1.58 x 10 ⁴ TCID ₅₀ /ml
	Measles	1.58 x 10 ⁴ TCID ₅₀ /ml
	Influenza B	3.16 x 10 ⁶ TCID ₅₀ /ml
	Influenza A H3N2	1 x 10 ⁵ TCID ₅₀ /ml
	Influenza A H1N1	3.16 x 10 ⁵ TCID ₅₀ /ml
Virus strain	Human Rhinovirus 2	2.81 x 10 ⁴ TCID5 ₅₀ /ml
virus strain	Human Rhinovirus 16	8.89 x 10 ⁶ TCID ₅₀ /ml
	Human Rhinovirus 14	1.58 x 10 ⁶ TCID ₅₀ /ml
	Human coronavirus OC43	1 x 10 ⁶ TCID ₅₀ /ml
	Human coronavirus NL63	1x 10 ⁶ TCID ₅₀ /ml
	Human coronavirus HKU1	1x 10 ⁶ TCID ₅₀ /ml
	Human coronavirus 229E	5x 10 ⁵ TCID ₅₀ /ml
	MERS coronavirus Florida	1.17x10 ⁴ TCID ₅₀ /ml
	Adenovirus type 7	1.58 x 10 ⁵ TCID ₅₀ /ml
	Adenovirus type 3	3.16 x 10 ⁴ TCID ₅₀ /ml
	Arcanobacterium	1.0x10 ⁸ org/ml
	Corynebacterium	1.0x10 ⁸ org/ml
	Escherichia coli	1.0x10 ⁸ org/ml
	Moraxella catarrhalis	1.0x10 ⁸ org/ml
	Neisseria lactamica	1.0x10 ⁸ org/ml
	Neisseria subflava	1.0x 10 ⁸ org/ml
Bacteria	Pseudomonas aeruginosa	1.0x10 ⁸ org/ml
	Staphylococcus aureus subspaureus	1.0x10 ⁸ org/ml
	Staphylococcus epidermidis	1.0x10 ⁸ org/ml
	Streptococcus pneumoniae	1.0x10 ⁸ org/ml
	Streptococcus pygenes	1.0x10 ⁸ org/ml
	Streptococcus salivarius	1.0x10 ⁸ org/ml
	Streptococcus sp group F	1.0x10 ⁸ org/ml
Yeast	Candida albicans	1.0x10 ⁸ org/ml

4. Interference Substances

The following substances that were evaluated had no impact on ACROTM SARS-CoV-2 Antigen Rapid Test (Swab) at concentrations listed below.

Substance	Concentration	
Whole Blood	20μl/ml	
Mucin	50μg/ml	
Budesonide Nasal Spray	200μl/ml	
Dexamethasone	0.8mg/ml	
Flunisolide	6.8ng/ml	
Mupirocin	12mg/ml	
Oxymetazoline	0.6mg/ml	
Phenylephrine	12mg/ml	
Rebetol	4.5μg/ml	
Relenza	282ng/ml	
Tamiflu	1.1µg/ml	
Tobramycin	2.43mg/ml	

5. Repeatability and Reproducibility

Repeatability and Reproducibility were evaluated using 3 standard controls. There were no differences observed within-run, between-run, between-lots, between-sites and between days.

1. Westgard JO, Barry PL, Hunt MR, Groth T. A multi-rule Shewhart for quality control in clinical chemistry, Clinical Chemistry 1981;27:493-501

Index of Symbols

		muex o	Symbols
IVD	For in vitro diagnostic use only	Σ	Tests per kit
2°C - 30°C	Store between 2-30°C	\square	Use by
®	Do not use if package is damaged	LOT	Lot Number
~	Manufacturer		Consult Instructions For Use

	EC REP	Authorized		
١		Representative		
	2	Do not reuse		
	REF	Catalog #		
or				



ACRO BIOTECH, Inc.

9500 Seventh Street, Unit M, Rancho Cucamonga, CA 91730, U.S.A.



Sterile swabs:



Jiangsu Changfeng Medical Industry Co., Ltd Touqiao Town, Guangling District, Yangzhou,Jiangsu 225109 China



Copan Italia S.p.A. Via F. Perotti, 10 25125 Brescia -Italy www.copangroup.com

C € 0123



Company LLC 31 School Street Guilford, Maine 04443-0149 207-876-3311 puritanmedproducts.com

€ 0086



Medico Technology Co., Ltd Room 201 of Building 14th and Buildin 17th Hengyi Lane, Yuanhu Road, Zhangbei industrial Park, Longcheng Street, Longgang district, Shenzhen, Guangdong, China www.medicoswab.com

C € 0413

EC REP

Number: 146370305 Effective Date: 2021-04-06