

PRODUCT CATALOGUE- FLUORESCENCE IMMUNOASSAY



Fluorescence Immunoassay Fluorescence Immunoassay Analyzers

Fluorescence Immunoassay Analyzer (AFR-ACO600)

Technical Specifications

Methodology :	Fluorescence Immunoassay
Test Format :	Cassette
Measurement :	Quantitative, Qualitative
Reading Time :	<20 Seconds
Display :	7" LCD touch color screen
Memory :	10000 Records
Dimensions(L*W*H) :	233*238*229 mm
Weight :	2.3 Kg
Storage Conditions :	5~45 °C, RH ≤85%
Working Conditions :	10~30 °C, RH ≤85%
Printer :	Built-in thermal printer
Power Source :	100-240 VAC,50/60Hz
Compatibility :	LIS
Scanner :	Laser barcode scanner with USB port
Connectivity :	USB 2 ports, LAN port, RS232 COM port (WLAN is optional functions)



Fluorescence Immunoassay Analyzer (ACA201)

Technical Specifications

Methodology	Fluorescence Immunoassay
Test Formats	Cassette
Measurement	Quantitative, Qualitative
Read Time	< 20 seconds
Dimensions (L*W*H)	244.5mm*240mm*201.4mm
Display Dimensions	7 inches
Weight	2.2 kg
OS	Windows-Compatible proprietary program
Power Supply	DC 24 V 4A
Adapter Input	AC 100~240 V, 50/60 Hz, 60 W
Connecting Ports	COM, USB, LAN
Memory	10000 records
Operating Conditions	10-30 °C
Relative Humidity	RH ≤ 85%
Atmospheric pressure	86 kPa~106 kPa
Printer	Built-in Thermal printer



Fluorescence Immunoassay Analyzer - 12 Channels (ACA202)

Technical Specifications

Methodology	Fluorescence Immunoassay
Test Formats	Cassette
Measurement	Quantitative, Qualitative
Read Time	< 20 seconds
Dimensions (L*W*H)	405mm*500mm*265mm
Display Dimensions	10.1 inch color LCD
Weight	15 kg
Maximum Input Power	200VA
Working Voltage	AC 100-240V 50/60 Hz
Connecting Ports	COM*1 USB *4 LAN * 1
Memory	10000 records
Operating Conditions	10-30 °C
Relative Humidity	RH ≤ 85%
Excitation / Emission band	365/610 (nm)
Printer	Built-in Thermal printer
Test channels	12



Fluorescence Immunoassay Portable Analyzer (ACA203)

Technical Specifications

Methodology :	Fluorescence Immunoassay
Test formats :	Cassette
Measurement :	Quantitative, Qualitative
Read Time :	< 20 seconds
Dimension :	275mm*105mm*100mm (L*W*H)
Screen size :	5 inches
Weight :	1kg
Spectrum :	Excitation spectrum: Mean wavelength $\lambda_0=365\text{nm}$ Reception spectrum: Mean wavelength $\lambda_1=610\text{nm}$
Memory :	10000 records
Printer :	Built-in thermal printer
Ports :	LAN×1, USB×1, Type-C×1



Fluorescence Immunoassay

Fluorescence Immunoassay Tests

Cat. No.	Product	Specimen	Cut-off	Test Range	Pack	CE status
Inflammation Marker						
ACO-CRP-402	CRP	Finger prick WB, WB/S/P	hsCRP:< 1 mg/L, CRP:<10 mg/L	0.5-200 mg/L	10T/25T	CE
ACO-MXCR-402	MXA/CRP	Finger prick WB, WB/S/P	MxA:<15 ng/mL; CRP:<10 mg/L	MxA:10-300 ng/mL; CRP:1-200 mg/L	10T/25T	CE
ACO-PCT-402	PCT	Finger prick WB, WB/S/P	<0.5 ng/mL	0.1-50 ng/ml	10T/25T	CE
ACO-IL6-402	IL-6	Finger prick WB, WB/S/P	<10 pg/mL	3.0-5000 pg/ml	10T/25T	CE
ACO-SAA-402	SAA	Finger prick WB, WB/S/P	<10 mg/L	1-300 mg/L	10T/25T	CE
Tumor Marker						
ACO-PSA-302	PSA	S/P	<4 ng/mL	2-100 ng/mL	10T/25T	Non-CE
ACO-PSA-402	PSA	Finger prick WB, WB/S/P	<4 ng/mL	0.1-100 ng/mL	10T/25T	Non-CE
ACO-FPSA-302	FPSA	S/P	/	0.1-30 ng/mL	10T/25T	Non-CE
ACO-AFP-302	AFP	S/P	<20 ng/mL	5-400 ng/mL	10T/25T	CE
ACO-CEA-302	CEA	S/P	<4 ng/mL	1-500 ng/mL	10T/25T	CE
ACO-FOB-602	FOB	Feces	<50 ng/ml	10-1000 ng/ml	10T/25T	CE
ACO-CA211-402	CYFRA21-1	Finger prick WB, WB/S/P	<3.3 ng/mL	1.0-300ng/mL	10T/25T	CE
ACO-CA153-402	CA153	Finger prick WB, WB/S/P	< 30.0 U/mL	2.0-300.0 U/mL	10T/25T	CE
ACO-CA199-402	CA19-9	Finger prick WB, WB/S/P	<37 U/mL	2-2000 U/mL	10T/25T	CE
ACO-CA125-402	CA125	Finger prick WB, WB/S/P	<35 U/mL	5-600 U/mL	10T/25T	CE
Cardiac Marker						
ACO-CTI-402	cTnI	Finger prick WB, WB/S/P	<0.5 ng/mL	0.1-40 ng/mL	10T/25T	CE
ACO-SCTI-402	High-Sensitivity Cardiac Troponin I	WB/S/P	<0.3 ng/mL	0.02-20 ng/mL	10T/25T	CE
ACO-CTNT-402	cTnT	Finger prick WB, WB/S/P	<0.5 ng/ml	0.2-40 ng/ml	10T/25T	CE
ACO-SCTT-402	High-Sensitivity Cardiac Troponin T	Finger prick WB, WB/S/P	<0.1 ng/mL	0.01-10 ng/mL	10T/25T	CE
ACO-NBNP-402	NT-proBNP	Finger prick WB, WB/S/P	/	300-22000 pg/mL	10T/25T	CE
ACO-SNBNP-402	High-Sensitivity NT-proBNP	Finger prick WB, WB/S/P	/	30-30000 pg/mL	10T/25T	CE
ACO-CKMB-402	CK-MB	Finger prick WB, WB/S/P	<5 ng/mL	0.2-75 ng/mL	10T/25T	CE
ACO-HFA-402	H-FABP	Finger prick WB, WB/S/P	<8 ng/mL	1-120 ng/mL	10T/25T	CE
ACO-MYO-402	Myoglobin	Finger prick WB, WB/S/P	< 90 ng/mL	5-200 ng/mL	10T/25T	CE
ACO-AMI-435	Troponin I/ Myoglobin/ CK-MB(3 in 1)	Finger prick WB, WB/S/P	cTnI: <0.5 ng/mL MYO: <90 ng/mL CK-MB: <5 ng/mL	0.1-40 ng/mL 5-200 ng/mL 0.2-75 ng/mL	10T/25T	CE
ACO-AMIT-435	Troponin T/ Myoglobin/ CK-MB (3in1)	Finger prick WB, WB/S/P	cTnT: <0.2 ng/mL MYO: <90 ng/mL CK-MB: <5 ng/mL	0.1-40 ng/mL 5-200 ng/mL 0.2-75 ng/mL	10T/25T	CE
ACO-CPC-435	Troponin I/ D-Dimer/ NT-proBNP (3 in 1)	Finger prick WB, WB/S/P	cTnI: <0.5 ng/mL D-Dimer: <500 ng/mL NT-proBNP: <0.45 ng/mL	0.1-40 ng/mL 100-10000 ng/mL 0.3-22 ng/mL	10T/25T	CE
ACO-LP-402	Lp-PLA2	Finger prick WB, WB/S/P	250ng/mL	10-1000 ng/mL	10T/25T	CE
ACO-ST2-402	ST2	Finger prick WB, WB/S/P	<35 ng/mL	1-500 ng/mL	10T/25T	CE
ACO-MPO-402	MPO	Finger prick WB, WB/S/P	350 ng/mL	30-1500 ng/mL	10T/25T	CE
Coagulation						
ACO-DDM-402	D-Dimer	Finger prick WB, WB/P	<0.5 mg/L	0.1-10 mg/L	10T/25T	CE
acrFDP402F	FDP	Finger prick WB, WB/S/P	<5 mg/L	2-80 mg/L	10T/25T	CE
Endocrinology						
ACO-BHCG-402	-hCG	Finger prick WB, WB/S/P	< 5 mIU/mL	2-200,000 mIU/mL	10T/25T	CE
ACO-FSH-402	FSH	Finger prick WB, WB/S/P	/	1-100 mIU/mL	10T/25T	CE
ACO-LH-402	LH	Finger prick WB, WB/S/P	/	1-300 mIU/mL	10T/25T	CE
ACO-T3-402	T3	Finger prick WB, WB/S/P	1.23-3.08 nmol/L	0.62-9.24 nmol/L	10T/25T	CE
ACO-FT3-402	FT3	Finger prick WB, WB/S/P	3.1-6.8 pmol/L	1.5-46 pmol/L	10T/25T	CE
ACO-T4-402	T4	Finger prick WB, WB/S/P	65-155 nmol/L	12.87-300 nmol/L	10T/25T	CE
ACO-FT4-402	FT4	Finger prick WB, WB/S/P	9-19.05 pmol/L	1-100 pmol/L	10T/25T	CE
ACO-TES-302	Testosterone	S/P	FeMale:0-2 ng/mL, Male:2.6-12 ng/mL	1.0-20 ng/mL	10T/25T	CE
ACO-TSH-402	TSH	Finger prick WB, WB/S/P	0.27-4.2 µIU/mL	0.1-200 µIU/mL	10T/25T	CE
ACO-P4-302	Progesterone (P4)	S/P	/	1.4-60 ng/mL	10T/25T	CE
acrP4402F	Progesterone (P4)	Finger prick WB, WB/S/P	/	1-60 ng/mL	10T/25T	CE
ACO-FER-402	Ferritin	Finger prick WB, WB/S/P	Male:30-400 ng/mL, Female:13-150 ng/mL	1-2000 ng/mL	10T/25T	CE
ACO-AMH-402	AMH	Finger prick WB, WB/S/P	/	0.05 - 25 ng/ml	10T/25T	CE
ACO-PRL-402	PRL	Finger prick WB, WB/S/P	Female (Non-pregnant): 2- 25 ng/mL Female (Pregnant): 10 - 209 ng/mL Male:2 - 18 ng/mL	1.0 - 200 ng/mL	10T/25T	CE
ACO-COR-302	Cortisol	S/P	140 - 690 nmol/L	50 nmol/L - 1000 nmol/L	10T/25T	CE
ACO-BHCG-202	Beta-HCG	Finger prick WB, WB/S/P/Urine	S/P: <5 mIU/mL,Urine :<25 mIU/mL	5-100,000 mIU/mL	10T/25T	CE
ACO-E2-402	Estradiol (E2)	Finger prick WB, WB/S/P	/	10-3000 pg/mL	10T/25T	CE
ACO-TG-402	Thyroglobulin (TG)	Finger prick WB, WB/S/P	3.5-77 ng/mL	1.0-500 ng/mL	10T/25T	CE
ACO-TGB-402	TGAb	Finger prick WB, WB/S/P	115 IU/mL	10-4000 IU/mL	10T/25T	CE
ACO-PTH-402	PTH	Finger prick WB, WB/S/P	0.015 ng/mL-0.065 ng/mL	0.01-5 ng/mL	10T/25T	CE
ACO-TPOB-402	TPOAb	Finger prick WB, WB/S/P	34 IU/mL	5-600 IU/mL	10T/25T	CE
ACO-INS-402	Insulin	Finger prick WB, WB/S/P	2.6-24.9 mU/L	1.0-1000 mU/L	10T/25T	CE
ACO-IGF-502	iGFBP-1	Finger prick WB, WB/S/P	<20 ng/mL	5-1000 ng/mL	10T/25T	CE
ACO-FFN-502	Fetal Fibronectin (fFN)	Vaginal Secretion	<50 ng/mL	5-1000 ng/mL	10T/25T	CE
acrGFN502F	iGFBP-1/fFN Combo	Vaginal secretion	iGFBP-1:<20 ng/mL; fFN:<50 ng/mL	5-1000 ng/mL	10T/20T	CE
ACO-CP-402	C-peptide	Finger prick WB, WB/S/P/Urine	WB/S/P: 0.8-4.2 ng/mL, Urine:36.4-81.36 µg/24h	WB/S/P: 0.04-40.0 ng/mL, Urine:0.16-160.0 ng/mL	10T/25T	CE
ACO-HGH-402	hGH	Finger prick WB, WB/S/P	/	0.1-50.0 ng/mL	10T/25T	CE

Cat. No.	Product	Specimen	Cut-off	Test Range	Pack	CE status
Nephrology						
ACO-ALB-102	mAlb	Urine	< 20 mg/L	5-300 mg/L	10T/25T	CE
ACO-CYSC-402	CysC	Finger prick WB, WB/S/P	<1.1 mg/L	0.1-9.0 mg/L	10T/25T	CE
ACO-BMG-402	2MG	Finger prick WB, WB/S/P	<2.2 mg/L	0.2-20 mg/L	10T/25T	CE
ACO-NGAL-102	N-GAL	Urine	<10 ng/mL	10-1500 ng/mL	10T/25T	CE
ACO-CRE-402	Creatinine	Finger prick WB, WB/S/P	Female: 45-106 µmol/L; Male: 59-115 µmol/L	5-2000 µmol/L	10T/25T	CE
ACO-ACR-102	Microalbumin/ Creatinine Ratio	Urine	<2.5 mg/mmol	0.1-300 mg/mmol	10T/25T	CE
Diabetes						
ACO-HBA-402	HbA1c	WB	4-6%	4-14.5%	10T/25T	CE
Infectious Diseases						
ACO-STA-502	Strep A	Swab	<1 COI	Qualitative (0-100 COI)	10T/20T	CE
ACO-FLU-502	FLU A+B	Swab	<1 COI	Qualitative (0-100 COI)	10T/20T	CE
ACO-RSV-502	RSV	Swab	<1 COI	Qualitative (0-100 COI)	10T/20T	CE
ACO-NCP-502	COVID-19 Antigen	Swab	<1 COI	Qualitative (1-1000 COI)	20T	CE
ACO-NCP-402	COVID-19 IgG/IgM	Finger prick WB, WB/S/P	<1 COI	Qualitative (0-20 COI)	10T/25T	CE
ACO-CSG-402	COVID-19 S-RBD IgG Antibody	Finger prick WB, WB/S/P	<1 COI	Semi-quantitative (0.5-100 COI)	10T/25T	CE
ACO-CIC-525	COVID-19 Ag/Influenza A+B Combo	NP Swab	<1 COI	Qualitative (COVID-19 Ag: 0-1000 COI FLUA/FLUB: 0-100 COI)	20T	CE
ACO-DEN-402	Dengue IgG/IgM	Finger prick WB, WB/S/P	<1 COI	Qualitative (0-100 COI)	10T/25T	CE
ACO-DES-402	Dengue NS1	Finger prick WB, WB/S/P	<1 COI	Qualitative (0-100 COI)	10T/25T	CE
ACO-MAL-435	Malaria P.f./P.v./Pan	Finger prick WB, WB	<1 COI	Qualitative (0-20 COI)	10T/25T	CE
ACO-HCV-402	HCV	Finger prick WB, WB/S/P	<1 COI	Qualitative (0-100 COI)	10T/25T	Non-CE
ACO-HBSB-402	HBsAb	WB/S/P	< 10mIU/mL	5-1000 mIU/mL	10T/25T	Non-CE
ACO-HBSG-402	HBsAg	WB/S/P	<0.5 IU/mL	0.5-200 IU/mL	10T/25T	Non-CE
acrHBEB402F	HBeAb	Finger prick WB, WB/S/P	> 1.00 COI	Qualitative (0-100 COI)	10T/25T	Non-CE
acrHBEG402F	HBeAg	Finger prick WB, WB/S/P	< 1.00 COI	Qualitative (0-100 COI)	10T/25T	Non-CE
acrHBCB402F	HBeAb	Finger prick WB, WB/S/P	> 1.00 COI	Qualitative (0-100 COI)	10T/25T	Non-CE
ACO-HIV-402	HIV1/2	Finger prick WB, WB/S/P	<1 COI	Qualitative (1-20 COI)	10T/25T	Non-CE
ACO-HIG-402	HIV P24 Antigen	Finger prick WB, WB/S/P	<1 COI	Qualitative (0-100 COI)	10T/25T	Non-CE
ACO-SYP-402	Syphilis	Finger prick WB, WB/S/P	<1 COI	Qualitative (0-100 COI)	10T/25T	CE
ACO-TET-402	Tetanus IgG	Finger prick WB, WB/S/P	<0.1 ng/mL	0.05-10IU/mL	10T/25T	CE
ACO-HP-602	H. pylori Antigen	Feces	< 2.5 ng/mL	2-1500 ng/mL	10T/25T	CE
ACO-ROT-602	Rotavirus	Feces	<1 COI	Qualitative (0-100 COI)	10T/25T	CE
ACO-HSV1-402	HSV 1 IgG/IgM	Finger prick WB, WB/S/P	<1 COI	Qualitative (0-100 COI)	10T/25T	CE
ACO-HSV2-402	HSV 2 IgG/IgM	Finger prick WB, WB/S/P	<1 COI	Qualitative (0-100 COI)	10T/25T	CE
acrCMV402F	CMV	Finger prick WB, WB/S/P	<1 COI	Qualitative (0-100 COI)	10T/25T	Non-CE
ACO-RUB-402	Rubella IgG/IgM	Finger prick WB, WB/S/P	<1 COI	Qualitative (0-100 COI)	10T/25T	Non-CE
ACO-TOX-402	Toxo IgG/IgM	Finger prick WB, WB/S/P	<1 COI	Qualitative (0-100 COI)	10T/25T	Non-CE
ACO-MP-502	Mycoplasma pneumoniae Antigen	Throat Swab	<1 COI	Qualitative (0-100 COI)	10T/20T	CE
ACO-SP-102	Streptococcus pneumoniae	Urine	<1 COI	Qualitative (0-100 COI)	10T/25T	CE
ACO-LEG-102	Legionella pneumophila	Urine	<1 COI	Qualitative (0-100 COI)	10T/25T	CE
ACO-AD-502	Adenovirus pneumoniae Antigen	Nasopharyngeal Swab	<1 COI	Qualitative (0-100 COI)	10T/20T	CE
ACO-AR-502	Adenovirus pneumoniae and RSV Combo	Nasopharyngeal Swab	<1 COI	Qualitative (0-100 COI)	10T/20T	CE
ACO-HMPV-502	HMPV	Nasopharyngeal Swab	<1 COI	Qualitative (0-100 COI)	10T/20T	CE
ACO-HPIV-502	HPIV	Nasopharyngeal Swab	<1 COI	Qualitative (0-100 COI)	10T/20T	CE
ACO-TYP-402	Typhoid	Finger prick WB, WB/S/P	<1 COI	Qualitative (0-100 COI)	10T/25T	CE
ACO-NOV-602	Norovirus	Feces	<1 COI	Qualitative (0-100 COI)	10T/25T	CE
ACO-AD-602	Adenovirus Antigen	Feces	<1 COI	Qualitative (0-100 COI)	10T/25T	CE
ACO-RAC-602	Rotavirus and Adenovirus	Feces	<1 COI	Qualitative (0-100 COI)	10T/25T	CE
ACO-CDG-602	Clostridium difficile GDH	Feces	<1 COI	Qualitative (0-100 COI)	10T/25T	CE
ACO-CDT-625	Clostridium difficile Toxin A/Toxin B Combo	Feces	<1 COI	Qualitative (0-100 COI)	10T/25T	CE
ACO-CD-602	Clostridium difficile GDH+Toxin A/B	Feces	<1 COI	Qualitative (0-100 COI)	10T/25T	CE
ACO-MXG-402	Monkeypox Virus Antigen Filariasis IgG/IgM	Finger prick WB, WB/S/P /Rash Exudate Swab /Throat Swab	<1 COI	Qualitative (0-100 COI)	20T	CE
ACO-FIL-402	Leptospira IgG/IgM	Finger prick WB, WB/S/P	<1 COI	Qualitative (0-100 COI)	10T/25T	CE
ACO-LEP-402	Leptospira IgG/IgM	Finger prick WB, WB/S/P	<1 COI	Qualitative (0-100 COI)	10T/25T	CE
ACO-TSU-402	Scrub Typhus IgG/IgM	Finger prick WB, WB/S/P	<1 COI	Qualitative (0-100 COI)	10T/25T	CE
ACO-ASO-402	Anti-streptolysin O (ASO)	Finger prick WB, WB/S/P	Adults: 200 IU/mL Children: <250 IU/mL	25-800 IU/mL	10T/25T	CE
Autoimmunity						
ACO-IGE-402	IgE	Finger prick WB, WB/S/P	100 IU/ml	20-1600 IU/mL	10T/25T	CE
ACO-IGA-402	IgA	Finger prick WB, WB/S/P	0.7-4.3 g/L	0.1-7.0 g/L	10T/25T	CE
ACO-RFM-402	Rheumatoid Factor IgM	Finger prick WB, WB/S/P	<20 IU/mL	8-200 IU/mL	10T/25T	CE
ACO-CCP-402	Anti-CCP	Finger prick WB, WB/S/P	<17.0 U/mL	3.0 U/mL- 300.0 U/mL	10T/25T	CE
Drugs of Abuse						
ACO-AMP-102	Single Drug - Amphetamine (AMP)	Urine	<1000 ng/mL	300-3000 ng/mL	10T/25T	CE
ACO-BAR-102	Single Drug - Barbiturates (BAR)	Urine	<300 ng/mL	50-1200 ng/mL	10T/25T	CE
ACO-BZO-102	Single Drug - Benzodiazepines (BZO)	Urine	<300 ng/mL	150-1500 ng/mL	10T/25T	CE
ACO-COC-102	Single Drug - Cocaine (COC)	Urine	<300 ng/mL	50-3000 ng/mL	10T/25T	CE
ACO-MET-102	Single Drug - Methamphetamine (MET)	Urine	<1000 ng/mL	100-5000 ng/mL	10T/25T	CE
ACO-MTD-102	Single Drug - Methadone (MTD)	Urine	<300 ng/mL	150-1500 ng/mL	10T/25T	CE
ACO-OPI-102	Single Drug - Opiate (OPI)	Urine	<2000 ng/mL	300-6000 ng/mL	10T/25T	CE
ACO-THC-102	Single Drug - Marijuana (THC)	Urine	<50 ng/mL	25-150 ng/mL	10T/25T	CE
ACO-APVP-102	Single Drug - -Pyrrolidinovalerophenone (-PVP)	Urine	<300 ng/mL	150-3000 ng/mL	10T/25T	CE
ACO-MDMA-102	Single Drug - Methylenedioxy-methamphetamine (MDMA)	Urine	<500 ng/mL	250-1500 ng/mL	10T/25T	CE
ACO-MOP-102	Single Drug - Morphine (MOP)	Urine	<300 ng/mL	150-6000 ng/mL	10T/25T	CE
ACO-ETO-H902	Single Drug - Etomidate (ETO)	Hair	0.2 ng/mg	Qualitative (0.2-1 ng/mg)	20T	CE
ACO-KET-H902	Single Drug - Ketamine (KET)	Hair	0.2 ng/mg	Qualitative (0.2-1 ng/mg)	20T	CE
ACO-MET-H902	Single Drug - Methamphetamine (MET)	Hair	0.2 ng/mg	Qualitative (0.2-1 ng/mg)	20T	CE
ACO-MOP-H902	Single Drug - Morphine (MOP)	Hair	0.2 ng/mg	Qualitative (0.2-1 ng/mg)	20T	CE
ACO-TML-H902	Single Drug - Tramadol (TML)	Hair	0.2 ng/mg	Qualitative (0.2-1 ng/mg)	20T	CE
Others						
ACO-VD-402	Vitamin D	Finger prick WB, WB/S/P	>30 ng/mL	5-100 ng/mL	10T/25T	CE
ACO-CAL-602	Calprotectin	Feces	<50 µg/g	1-1000 µg/g	10T/25T	CE
ACO-VB12-302	Vitamin B12	S/P	197-771 pg/mL or 145-569 pmol/L	50-2000 pg/mL or 36.9-1476 pmo/L	10T/25T	CE
ACO-FA-302	Folate	S/P	>3ng/mL	1-20 ng/mL	10T/25T	CE

ABOUT US

Acro Biotech is a California-based manufacturer of rapid tests offering the products from state-of-the-art ISO13485:2016 accredited CMDCAS-compliant manufacturing facilities. Our core team comprises of seasoned professionals in the field of rapid tests with more than two decades of experience in reputed rapid test manufacturing companies, working in senior positions. Our product range includes USFDA (510)K approved DOA tests as well as an expansive range comprising of women's health related rapid tests, cardiac markers, tumor markers, infectious disease and veterinary tests for international market.

Our tests work with a wide range of specimens, including urine, feces, oral fluid, whole blood, serum, plasma, hair and swabs.

At the cusp of innovation and commercialization, we are working hard to introduce new products in the market that would set the future growth path, including the combination of digitization and ease of tests with lateral flow readers and fluorescence immunoassay analyzers.

